

# SCIENTIFIC REPORT

## Grandio – Recommended polishers

VOCO GmbH, Knowledge Communication Department

Anton-Flettner-Str. 1-3  
27472 Cuxhaven, Germany

Tel.: +49 (0)4721-719-1111  
Fax: +49 (0)4721-719-109

info@voco.de  
www.voco.com



A hard material such as Grandio places special demands on the polisher. VOCO, therefore, tested all popular polishing systems and compiled a list of the most effective and fastest products for polishing Grandio. In addition to "Easy Gloss" from VOCO, the following systems are suitable for achieving an excellent polishing result in a short amount of time.

Name	Manufacturer	Manufacturer's Article No.
<b>Easygloss</b>	<b>VOCO</b>	<b>2050 (cup, 6 units)</b>
Ceramaster	Shofu	0121 (bullet), 0122 (lenticular), 0123 (cup) 0124 (mini-point)
DiaGloss	Edenta	19032RA (large point), 19033RA (small point) 19034RA (large cup), 19035RA (medium cup) 19038RA (disc), 19039RA (small cup)
Diapol	EVE Ernst Vetter	W11D (small point), W16D (small flame), W2D (large flame), W17D (small cup) W7D (large cup), W18D (lenticular) W17D (disc)
Identoflex Porzellan/Keramik NG	Kerr HAWE	ID 7251/12 (disc), ID 7261/12 (lenticular) ID 7271/12 (flame), ID 7281/12 (cup) ID 7291/12 (point)
Luster	Hager & Meisinger	2660 (assortment, 9 units) 9780 (flame), 9781 (cup) 9782 (disc), HM135 (cone) HM379 (flame), HM246 (point)
Occlubrush	Kerr HAWE	2520 (assortment, 15 units) 2503 (cup), 2504 (small cup), 2505 (point)
OptiShine	Kerr HAWE	2513 (cup, 3 units), 2514 (cup, 10 units)
PoGo	Dentsply	662025Y (assortment, 30 units)

All products can be sterilised (Last updated 12/2008).