



Traypurol®

PROVEN CLEANING CONCENTRATE



Traypurol®

PROVEN CLEANING CONCENTRATE

What is Traypurol®?

Traypurol is the ideal aide for the daily cleaning of instruments and impression trays in your surgery! Traypurol reliably removes residues of alginates, glass ionomer, carboxylate, phosphate and zinc oxide eugenol cements and is suitable for all non-oxidising metals and plastics. Moreover, it is biodegradable and provides gentle cleaning action.

How is Traypurol® used?

Using Traypurol is quite simple: It is added to water and dependably cleans trays and instruments placed in the solution. The material residues are partially dissolved by the Traypurol solution and can then be easily rinsed off under running water or removed with a brush. Traypurol can also be used in ultrasonic devices.

In which forms of presentation is Traypurol® available?

Traypurol is available in the proven 1 litre bottle and in an economical 5 litre canister for increased demand.

Your advantages

- Faster-acting
- Improved dissolving power
- Pleasant, subtle scent

Presentation

REF 2289 Bottle 1 I
REF 2290 Canister 5 I

REF 2304 Dispenser pump – canister 5 I



VOCO GmbH Anton-Flettner-Straße 1-3 27472 Cuxhaven Germany

Freecall: 00 800 44 444 555 Fax: +49 (0) 4721-719-140

info@voco.com www.voco.dental

Available from:		

