

## VOCO's composites – Classification

VOCO GmbH, Knowledge Communication Department

Anton-Flettner-Str. 1-3  
27472 Cuxhaven, Germany

Tel.: +49 (0)4721-719-1111  
Fax: +49 (0)4721-719-109

info@voco.de  
www.voco.com

Product	Material class	Filler degree [% by weight]
Admira	Micro-hybrid ORMOCER®	78
Admira Flow	Micro-hybrid ORMOCER®	63
Admira Fusion	Nano-hybrid ORMOCER®	84
Admira Fusion Flow	Nano-hybrid ORMOCER®	74
Admira Fusion x-tra	Nano-hybrid ORMOCER®, 4 mm	84
Admira Fusion x-base	Nano-hybrid ORMOCER®, 4 mm	72
Amaris	Micro-hybrid	80
Amaris Flow	Micro-hybrid	64
Amaris Gingiva	Micro-hybrid	80
Arabesk	Micro-hybrid	77
Arabesk Flow	Micro-hybrid	64
Arabesk Top	Micro-hybrid	77
Glasiosite	Compomer	78
Grandio	Nano-hybrid	87
Grandio Flow	Nano-hybrid	80
GrandioSO	Nano-hybrid	89
GrandioSO Heavy Flow	Nano-hybrid	83
GrandioSO Flow	Nano-hybrid	81
GrandioSO Light Flow	Nano-hybrid	76
GrandioSO x-tra	Nano-hybrid, 4 mm	86
Twinky Star	Compomer	78
Twinky Star Flow	Compomer	65
VisCalor bulk	Nano-hybrid, 4 mm	83
x-tra base	Micro-hybrid, 4 mm	75
x-tra fil	Micro-hybrid, 4 mm	86