

Les matériaux de restauration de VOCO – classification

VOCO GmbH, Département Communication-Formation

Anton-Flettner-Str. 1-3

D-27472 Cuxhaven

Tel.: +49 (0)4721-719-1111

Fax: +49 (0)4721-719-109

info@voco.de

www.voco.dental

Produit	Classification	Taux de charges [% en masse]
Admira	ORMOCER® microhybride	78
Admira Flow	ORMOCER® microhybride	63
Admira Fusion	ORMOCER® nanohybride	84
Admira Fusion Flow	ORMOCER® nanohybride	74
Admira Fusion x-tra	ORMOCER® nanohybride , 4 mm	84
Admira Fusion x-base	ORMOCER® nanohybride , 4 mm	72
Amaris	Microhybride	80
Amaris Flow	Microhybride	64
Amaris Gingiva	Microhybride	80
Arabesk	Microhybride	77
Arabesk Flow	Microhybride	64
Arabesk Top	Microhybride	77
Glasiosite	Compomère	78
Grandio	Nanohybride	87
Grandio Flow	Nanohybride	80
GrandioSO	Nanohybride	89
GrandioSO Heavy Flow	Nanohybride	83
GrandioSO Flow	Nanohybride	81
GrandioSO x-tra	Nanohybride, 4 mm	86
Twinky Star	Compomère	78
Twinky Star Flow	Compomère	65
x-tra base	Microhybride, 4 mm	75
x-tra fil	Microhybride, 4 mm	86