

VOCO HIGHLIGHTS

Practice

Valid from 01.04. - 30.09.2020

Highlight 1: VisCalor® bulk

Thermoviscous bulk fill composite

- Heating of the material makes it flowable for the application and then sculptable immediately afterwards (thermo-viscous-technology)
- Optimal flowing to margins and undercut regions – minimises risk of marginal gaps
- 4 mm bulk fill without covering layers
- Bubble-free application
- One universal and three aesthetic shades
- Slender cannula – ideal for hard-to-reach areas

Set VisCalor inclusive Dispenser

Caps 80 × 0.25 g (16 × universal, 16 × A1, 16 × A2, 32 × A3), VisCalor Dispenser

REF 6062

556,00 €



Preheating dispenser for composite Caps

VisCalor Dispenser – Preheating device

REF 9143

790,00 €

First
flowable,
then
sculptable

Caps 16 × 0.25 g universal
Caps 16 × 0.25 g A1
Caps 16 × 0.25 g A2
Caps 16 × 0.25 g A3

REF 6065
REF 6066
REF 6067
REF 6068

64,50 €

Highlight 2: Futurabond® U

Dual-curing universal adhesive

- One bond for all cases – makes all other adhesives in the practice redundant
- Self-etch, selective-etch or total-etch – free choice for the practitioner
- Outstanding versatility of application
- Moisture tolerant

SingleDose 50 pcs., accessories

REF 1571

97,50 €

SingleDose 200 pcs., accessories

REF 1572

341,00 €

Price
advantage
12%



Highlight 3: Futurabond® M+

Universal adhesive

- Can be applied with or without phosphoric acid etching
- Excellently suitable for direct or indirect restorations
- Secure adhesion to various materials such as metals, zirconium dioxide or aluminium oxide, as well as silicate ceramics, without additional primer

Bottle 5 ml, accessories

REF 1515

103,40 €

Bottle 3 × 5 ml

REF 1516

263,00 €

Price
advantage
15%



Highlight 4: VOCO Profluorid® Varnish

Fluoride-containing dental desensitising varnish (5 % NaF)

- Fast desensitisation and fluoride release (5 % NaF \triangleq 22,600 ppm fluoride)
- Easy to apply – moisture tolerant
- Available in six pleasant flavours
- Economical application in thin layers, high yield
- Esthetic, tooth-shaded varnish



SingleDose 50 x 0.40 ml

or

SingleDose 48 x 0.40 ml

(REF see table below)

SingleDose 200 x 0.40 ml

(REF see table below)

107,20 € ~~428,80 €~~ 316,00 €

save**
26%



| Presentations SingleDose | Melon | Caramel | Cherry | Mint | Bubble Gum | Cola Lime | Mixed*** |
|-----------------------------|-------|---------|--------|------|---------------|--------------|----------|
| 50 x 0.40 ml | 1269 | 2226 | 2227 | 2228 | 2238 | 2240 | - |
| 200 x 0.40 ml | 1271 | 2229 | 2230 | 2231 | 2239 | 2241 | - |
| 48 x 0.40 ml | - | - | - | - | - | - | 2225 |

***Mixed: each of 8 x Melon, Caramel, Cherry, Mint, Bubble Gum, Cola Lime

Highlight 5: Provicol® QM Aesthetic

Translucent temporary luting cement with calcium hydroxide, eugenol-free

- Highly aesthetic, recreates the appearance of enamel
- Increased strength for secure adhesion
- Pinpoint stability allows precise removal of excess material

QuickMix syringe 5 ml, accessories REF 1108

39,40 €

Set II QuickMix syringe 3 x 5 ml
(1 x Provicol QM, 1 x Provicol QM Plus,
1 x Provicol QM Aesthetic), accessories

REF 1109

75,00 €

NEW

save**
34%



Highlight 6: Ionolux®

Light-curing glass ionomer restorative material

- Aesthetic restorations without conditioning or bonding
- Excellent working time – individually adjustable by light-curing
- Does not stick to the instrument and can be modelled with ease

Application capsule 150 pcs. A3
+ 50 pcs. A3

REF 2989

499,00 €

save**
27%



Application capsule 50 pcs.
(10 x A1, 10 x A2, 30 x A3)
+ 20 pcs. A3

REF 2990

189,00 €

save**
29%



Highlight 7: Grandio®SO Light Flow

Low-viscosity nano-hybrid composite for delicate applications

- **Low viscosity** – outstanding flow properties and therefore optimally suitable for hard-to-reach areas and small cavities
- **Extra-fine cannula** – for precise and pinpoint application, dosage without excess residues
- **High performance** – excellent physical properties, e.g. filler content of 76% by weight
- **High aesthetics** – 8 shades for individual applications
- **Excellent handling** – NDT®-syringe – non-drip, compatible with all light-curing bonding materials

GrandioSO Light Flow – Set syringe 5 x 2 g
(A1, A2, A3, A3.5, WO), accessories

REF 2585

170,00 €
136,00 €

save**
20%



Syringe 2 x 2 g, accessories

| | | | |
|------|----------|---------|----------|
| A1 | REF 2586 | B1 | REF 2590 |
| A2 | REF 2587 | WO | REF 2591 |
| A3 | REF 2588 | OA2 | REF 2592 |
| A3.5 | REF 2589 | Incisal | REF 2593 |

68,00 €



Highlight 8: Admira® Fusion

Universal nano-hybrid ORMOCER® restorative material

- The worldwide first purely ceramic-based restorative material
- By far the lowest polymerisation shrinkage and particularly low level of shrinkage stress in comparison to all conventional restorative composites
- Inert, so excellent biocompatible and extremely resistant to discolouration
- Excellent handling, simple high-lustre polishing procedure coupled with high surface hardness guarantee first-class long-term results
- Compatible with all conventional bonding agents

The first pure ceramic for direct fillings.

Nothing but impressive advantages



Admira Fusion

Set + bond
syringe 5 × 3 g
(A2, A3, GA3.25, A3.5,
Admira Fusion x-tra),
shade guide + Futurabond U
SingleDose 20 pcs.

REF 2750 ~~336,90 €~~

save 30%

235,00 €



Set + bond
Caps 75 × 0.2 g
(15 each of A2, A3, GA3.25,
A3.5, Admira Fusion x-tra),
shade guide + Futurabond U
SingleDose 20 pcs.

REF 2780 ~~336,90 €~~

save 30%

235,00 €



4+1*

Syringe 3 g,
various shades

(REF see table below)

Caps 15 × 0.2 g,
various shades

62,70 €

Admira Fusion Flow

Set syringe 5 × 2 g
(A1, A2, A3, A3.5, WO),
shade guide, accessories

REF 2816 ~~203,90 €~~

save 30%

142,50 €



Syringe 2 × 2 g,
various shades

(REF see table below)

4+1*

74,20 €

THE ONLY REASON
NOT TO USE
ADMIRA FUSION
COULD BE YOU.

Admira Fusion (4 mm bulk fill) x-tra

Syringe 3 g universal

REF 2810

4+1*

Caps 15 × 0.2 g universal

REF 2811

62,70 €

Admira Fusion (flowable 4 mm bulk fill) x-base

Syringe 2 × 2 g universal,
accessories

REF 2812

4+1*

74,20 €

* If you order four packages via a VOCO dental consultant, you will receive another package for free. The free goods will be sent to you by VOCO.



| | A1 | A2 | A3 | A3.5 | A4 | B1 | B2 | B3 | C2 | D3 | OA1 | OA2 | OA3 | OA3.5 | BL | WO | Incisal | GA3.25 | GA5 | Mixed |
|------------------------------------|------|------|------|------|------|------|------|------|------|------|------|------|------|-------|------|------|---------|--------|------|-------|
| Admira Fusion Syringe 3 g | 2754 | 2755 | 2756 | 2758 | 2759 | 2762 | 2763 | 2764 | 2766 | 2768 | 2770 | 2771 | 2772 | 2773 | 2775 | - | 2776 | 2757 | 2760 | - |
| Admira Fusion Caps 15 × 0.2 g | 2782 | 2783 | 2784 | 2786 | 2787 | 2790 | 2791 | 2792 | 2794 | 2796 | 2798 | 2799 | 2800 | 2801 | 2803 | - | 2804 | 2785 | 2788 | 2806 |
| Admira Fusion Flow Syringe 2 × 2 g | 2818 | 2819 | 2820 | 2821 | 2822 | 2824 | 2825 | - | 2826 | - | - | 2829 | - | - | 2831 | 2832 | 2828 | - | - | - |

Highlight 9: V-Posil

Precision impression material, VPS

- Highly hydrophilic A-silicone for maximum precision
- Long working time combined with a short intraoral setting time
- Exceptional toughness and recovery after deformation ensure safety during or after removal
- Good hydrophilicity, even after setting, renders pouring and thus the prosthetic restoration perfect



Starter-Set Putty Fast / Light Fast

- | | |
|------------|----------------------------------|
| Putty Fast | Jar 2 x 450 ml, accessories |
| Light Fast | Cartridge 2 x 50 ml, accessories |

REF 2984

133,00 €

save**
15%



Starter-Set Putty Fast / Light Fast

- | | |
|------------|-----------------------------------|
| Putty Fast | Cartridge 2 x 380 ml, accessories |
| Light Fast | Cartridge 2 x 50 ml, accessories |

REF 2991

223,00 €

save**
15%



Starter-Set Heavy Soft Fast / Light Fast

- | | |
|-----------------|-----------------------------------|
| Heavy Soft Fast | Cartridge 2 x 380 ml, accessories |
| Light Fast | Cartridge 2 x 50 ml, accessories |

REF 2985

223,00 €

save**
15%



Starter-Set Mono Fast / Light Fast

- | | |
|------------|-----------------------------------|
| Mono Fast | Cartridge 2 x 380 ml, accessories |
| Light Fast | Cartridge 2 x 50 ml, accessories |

REF 2987

223,00 €

save**
15%



Highlight 10: Structur 3

Self-curing composite material for the fabrication of provisional crowns and bridges

- Quick, simple and reliable working, short setting time
- With natural fluorescence, available in 8 shades
- High compressive strength

Cartridge 50 ml, mixing tips type 6

- | | |
|---------------|-------------|
| A1 REF 2503 | B1 REF 2507 |
| A2 REF 2504 | B3 REF 2508 |
| A3 REF 2505 | C2 REF 2509 |
| A3.5 REF 2506 | BL REF 2510 |

106,00 €

Cartridge 5 x 50 ml, mixing tips type 6

- | |
|-------------|
| A2 REF 2512 |
| A3 REF 2513 |

~~530,00 €~~

485,00 €

save**
8%



Gloss without polishing



**compared with individual purchase

YOUR ORDER

Please contact your dental depot or your VOCO dental consultant to order the VOCO Highlights presented here. Offers are only valid at participating retail partners. Special prices offered by your depot are always taken into consideration. All prices are RRP and do not include VAT.

THIS BROCHURE WAS HANDED OUT BY: