

TOP PRODUCTS BY VOCO – RELIABLE AND INNOVATIVE



World's first composite with thermo-viscous-technology • Innovative material By pre-warming, the material is

Innovative material
 By pre-warming, the material is
 flowable during application and can be
 sculpted immediately afterwards
 (thermo-viscous-technology)

 Time-saving Efficient restorations with only one material

 Excellent handling Bubble-free application with a narrow and long cannula





... warm i



... flow it



... sculpt

VisCalor® bulk

REF 6062 Set Dispenser, Caps 80×0.25 g ($16 \times$ universal, $16 \times A1$, $16 \times A2$, $32 \times A3$), VisCalor Dispenser

REF 6063 Set Caps Warmer Caps 80×0.25 g (16×1.00 universal, 16×1.00 A1, 16×1.00 A2, 32×1.00 A3), Caps Warmer



Applying the heated VisCalor bulk with its superb flowable consistency and slender cannula

Adapting the sculptable VisCalor bulk



Polished and rehydrated class II cavities the day after filling placement

Source: Dr. Fernando Feitosa, Americana / Brazil

Celalux® 3

- State-of-the-art LED technology with high light intensity (approx. 1,300 mW / cm²)
- Wavelength range 450 480 nm
- Slim and ergonomic design (pen-shaped); weight of the handpiece, including device base and light guide, weighs only 70 g
- Indicates how much time has elapsed by vibrating at the start, and again after 10 s and 20 s, instead of the acoustic signals previously used
- Easy to handle, cordless
- Autoclavable 8 mm light guide
- The spare device base supplied ensures the device is always ready for use
- Battery is very easy to change

REF 9090 Celalux 3 (handpiece, battery charger, power cord with interchangeable mains adapters, 2 Device Bases, 8 mm light guide, rubber shield, operating instructions)



Take **FIVE**to match them all!

- All 16 VITA® classical shades can be covered with just 5 Cluster-Shades
- Only 10 seconds light-curing for all shades
- Excellent biocompatible, as no classic monomers are contained
- By far the lowest polymerisation shrinkage (1.25 vol.%)
- Compatible with all conventional bonding agents







Admira Fusion 5

VITA® classical

A1			A2			А3		A3.5				A4			
A1	В1	C1	D2	A2	B2	А3	В3	D4	A3.5	C2	С3	D3	B4	A4	C4

Admira® Fusion 5

REF 6243 Trial pack
10 × 0.2 g Caps (2 each of × A1,
A2, A3, A3.5, A4),
10 × Futurabond U SingleDose,
shade card, accessories
REF 6225 Set Syringe 5 × 3 g (A1, A2, A3,
A3.5, A4), shade card
REF 6235 Caps 75 × 0.2 g (15 each of ×
A1, A2, A3, A3.5, A4), shade
card



Admira® Fusion

REF 2750 Set + bond

syringe 5 × 3 g (A2, A3, GA3.25, A3.5, Admira Fusion x-tra), shade guide + Futurabond U SingleDose 20 pcs., accessories

REF 2780 Set + bond

Caps 75 × 0.2 g (15 each of A2, A3, GA3.25, A3.5, Admira Fusion x-tra), shade guide + Futurabond U SingleDose 20 pcs., accessories

Syringe 3 g¹ Caps 15×0.2 g¹

Admira® Fusion Flow

REF 2816 Set syringe 5 × 2 g (A1, A2, A3, A3.5, WO) shade guide, accessories

Syringe 2×2 g, various shades 1

¹ Further trading forms and single shades can be found under www.voco.dental or in the main catalogue

Admira® Fusion x-base

REF 2812 Syringe 2×2 g universal

Admira® Fusion x-tra

REF 2810 Syringe 3 g universal REF 2811 Caps 15×0.2 g universal




'I am convinced by GrandioSO's similarity to natural teeth!'

Grandio®SO

Our composite for high demands. The filling material scores with tooth-like properties combined with high user-friendliness and natural aesthetics.

- Universal nano-hybrid restorative material
- Excellent physical properties Very high filler content of 89 % w/w
- For the highest demands in both the anterior and posterior tooth regions



Grandio®SO Light Flow

- Flowable nano-hybrid restorative material - low viscous
- Precise application thanks to an extra-fine cannula and therefore optimally suitable for hard-to-reach areas and small cavities









- Flowable universal nano-hybrid restorative material - medium viscous
- Outstanding flow behaviour, universal application and precise application





Grandio®SO Heavy Flow

- Flowable universal nano-hybrid restorative material - highly viscous
- · Increased stability, i.e. no undesirable runoff and universal application







REF 2585 Set syringe 5×2 g (A1, A2, A3, A3.5, WO), shade guide, accessories

Syringe $2 \times 2 \text{ g}^1$

REF 2720 Set syringe 5×2 g (A1, A2, A3, A3.5, WO), shade guide, accessories

Syringe 2×2 g / Caps 16×0.25 g¹

REF 2680 Set syringe 5×2 g (A2, A3, A3.5, WO, GA5), shade guide, accessories

Syringe 2×2 g / Caps 16×0.25 g¹

¹Further trading forms and single shades can be found under www.voco.dental or in the main catalogue

Universal adhesive – It's your choice!

A bonding for all cases! With Futurabond U, the only dual-curing universal adhesive in *SingleDose*, we have a clear mission: To provide you a bonding which offers only advantages and reliability and makes any other adhesive system in your practice redundant.

Futurabond® U

- Dual-curing universal adhesive
- One bond for all cases makes all other adhesives in the practice redundant
- Self-etch, selective-etch or total-etch free choice for the practitioner
- Outstanding versatility of application



REF 1571 SingleDose 50 pcs., accessories

REF 1572 SingleDose 200 pcs., accessories

Self-etch



Selective-etch



Total-etch



Each drop a win! A universal adhesive that is right for every situation.

Futurabond® M+

- Universal adhesive
- Can be applied with or without phosphoric acid etching
- Excellently suitable for direct or indirect restorations
- Secure adhesion to various materials such as metals, zirconium dioxide or aluminium oxide, as well as silicate ceramics, without additional primer



REF 1515 Bottle 5 ml, accessories

REF 1516 Bottle 3 × 5 ml

Futurabond DC is the proven and popular self-etch-bond for the application under all light, self or dual-curing restorative materials – for a strong and permanent adhesion!

Futurabond® DC

- Dual-curing self-etch-bond reinforced with nano-fillers
- Chemical-curing in areas inaccessible with a curing light
- Futurabond DC achieves the same adhesive properties as total-etch preparations without separate etching of the tooth substance
- Moisture tolerant



REF 1164 SingleDose 50 pcs., accessories

REF 1165 SingleDose 200 pcs., accessories

Grandio - The classic nano-hybrid

With Grandio, VOCO has set standards in 2003 as the world's first supplier of a new generation of nano-hybrid composites. Proven millions of times on all continents, Grandio is still the nano-composite of choice! The innovative nano-technology makes Grandio a high-performance composite, which is known for its physical qualities as well as its optimal handling properties – everywhere.

Grandio®

- Ideal consistency and simplest handling
- E-modulus comparable to natural tooth substance
- Pronounced abrasion resistance and chewing stability
- High biocompatibilty thanks to low resin content



REF 1800 Set syringe 5×4 g (A1, A2, A3, A3.5, B2), shade guide

REF 1830 Set Caps 50 × 0.25 g (10 each of A1, A2, A3, A3.5, B2), shade guide

Further trading forms and single shades can be found under www.voco.dental or in the main catalogue



Initial situation
Source: Dr. Marcelo Balsamo



Preparation and lining with Grandio Flow



Finished restoration with Grandio


Grandio Flow – also unrivalled as flowable!

From the restoration of minimally invasive cavities to extended fissure sealing and the repair of veneers - with Grandio Flow, all the advantages of Grandio can also be used in the flowable version. Thanks to its clear technological advantage, Grandio Flow leaves conventional flowables far behind: The resin content is up to 50% lower than that of conventional flowables, which leads to top results in the areas of surface hardness and abrasion resistance.

Grandio® Flow


- Optimal flowability
- High filler content of 80 % w/w
- Exact shade match to Grandio
- Microhybrid composite like physical properties





Traditional flowable, unfilled area

REF 1860 Set syringe 8×2 g (A1, $2 \times A2$, $2 \times A3$, A3.5, B2, OA2), shade guide, accessories



Grandio Flow, filled with nano-particles

Further trading forms and single shades can be found under www.voco.dental or in the main catalogue

Amaris – for brilliant results

Naturally simple, simply beautiful! When it comes to the highest aesthetics, we recommend Amaris. The special thing about this composite of the extra class: The system allows you a high degree of flexibility in shading as the two coordinated working steps lets you change the shade continously – even during layering! Thus, a highly aesthetic result does not only depend on the initial shade selection.

Amaris®

- Optimal polishability, abrasion resistant and glossy
- Two easy steps, low material costs
- Nature-identical layering
- Covering of the entire shade range through the combination of only eight components



Further trading forms and single shades can be found under www.voco.dental or in the main catalogue

REF 1949 Starter set Caps 32×0.25 g (4 each of O1, O2, O3, O4, O5, TL, TN, TD), Amaris Flow syringe 2×1.8 g (HT, HO), shade guide, accessories

REF 1930 Set syringe 8×4 g (01, 02, 03, 04, 05, TL, TN, TD), Amaris Flow syringe 2×1.8 g (HT, HO), shade guide, accessories



Application of Amaris TL to reconstruct the palatal enamel portion

Source: Dr. Sânzio Marques – Passos, Brazil



Modelling with Amaris 01, natural appearance shaping of the mamelons with Amaris Flow HO



Completed anterior restoration

x-tra fil: Amalgam – no thanks!

You all might know the patients wish for an amalgam-free restoration but at amalgam prices. x-tra fil was specially designed for the basic treatment in the posterior region (class I and II) and for core build-ups. Thanks to the possibility of applying increment thicknesses of up to 4 mm in only one working step and the short light-curing time of 10 seconds, you and your patients benefit from a very short treatment time.

x-tra fil®

- Curing depth 4 mm
- Only 10 s curing time per layer
- Universal shade
- Reduced working time, especially in combination with any Futurabond SingleDose adhesive





REF 1739 Syringe 10×5 g

REF 1740 Syringe 5 g


REF 1741 Caps 20×0.25 g

The cord from the capsule

The VOCO Retraction Paste is your indispensable helper for the effective temporary widening and dry gingival sulcus. This creates the condition for successful conventional or digital impressions.

VOCO Retraction Paste

- Astringent paste for an effective temporary widening and dry gingival sulcus
- Thin cannula with flexible tip easy and pinpoint application into the sulcus
- Viscosity Change –
 paste consistency varies during
 application and sulcus widening
- Good visibility contrasty to the gingiva





REF 1012 Caps 25×0.3 g

REF 1013 Caps 100 × 0.3 g

Impressively precise!

Precisely fitted dentures require a material that ensures optimal impression taking with its properties – from the beginning to the pouring. This service concept has been consistently implemented in V-Posil.

V-Posil®

· Precision impression material, VPS

Highly hydrophilic A-silicone for maximum precision

 Long working time combined with a short intraoral setting time

 Exceptional toughness and recovery after deformation ensure safety during or after removal

 Good hydrophilicity, even after setting, renders pouring and thus the prosthetic restoration perfect



Starter-Set Putty Fast / Light Fast

REF 2991

Putty Fast Cartridge 2 \times 380 ml Light Fast Cartridge 2 \times 50 ml

Starter-Set Heavy Soft Fast / Light Fast

REF 2985

Heavy Soft Fast Cartridge 2×380 ml Light Fast Cartridge 2×50 ml

Starter-Set Mono Fast / Light Fast

REF 2987

Mono Fast Cartridge 2×380 ml Light Fast Cartridge 2×50 ml

Starter-Set Putty Fast / Light Fast

REF 2984

Putty Fast Jar $2 \times 450 \text{ ml}$ Light Fast Cartridge $2 \times 50 \text{ ml}$

Starter-Set Putty Soft Fast / Light Fast

REF 2904

Putty Soft Fast Jar 2×450 ml Light Fast Cartridge 2×50 ml

Linings for every day

A good lining is the base for a successful restoration. Whether for direct or indirect capping of the pulp or light-curing. With our linings you get to chose depending on the situation.

Ionoseal®

- Light-curing glass ionomer composite liner
- · Ready-to-use one-component material
- Time-saving: light-curing in seconds
- · Quick and hygienic application



REF 1326 Syringe 3×2.5 g, accessories



Prepared cavity
Source: Dr. Marcelo Balsamo



Direct application Ionoseal



Finished lining



Final restoration

Calcicur®

- Ready-to-use calcium hydroxide paste
- Suitable for use under any lining or restorative material
- · Immediately ready for use
- With high pH value (12.5) known for bacteriostatic effect



REF 1096 Syringe 3×2.5 g, accessories

Calcimol LC

- Light-curing radiopaque calcium hydroxide paste
- Effective protection of the pulp
- Ready-to-use one-component material
- Time-saving light-curing and direct application from the non-dripping NDT®-syringe



REF 1307 Syringe 2 x 2.5 g, accessories

REF 1047 Tube 2×5 g, accessories

Strong temporaries

Masticatory function, phonetics and aesthetics have to be right, even with a temporary restoration. With Structur, our temporary self-curing composite for the fabrication of high-quality provisional crowns and bridges, even durable longterm temporaries can be fabricated in just a few steps. The material is based on the known nano-hybrid technology from VOCO.

Structur 2

- Simple and fast working steps
- Less than 1 minute setting time in the mouth
- · Easy to shape and polish



REF 1492 Cartridge 5×75 g A2, mixing tips type 6 REF 1491 Cartridge 5×75 g A3, mixing tips type 6

Cartridge 75 g, accessories 1

Structur 3

- Processing
 - quick, simple and reliable thanks to the 1:1 cartridge system
 - short setting time
 - elastic phase allows simple removal from the prepared tooth
 - gloss without polishing only requires removal of the inhibition layer
- Durability
 - high compressive strength
 - excellent fracture resistance



REF 2512 Cartridge 5×50 ml A2, mixing tips type 6 REF 2513 Cartridge 5×50 ml A3,

mixing tips type 6

Cartridge 50 ml, accessories¹

¹Further trading forms and single shades can be found under www.voco.dental or in the main catalogue



Following removal during elastic phase: extraoral



Finished temporary bridge



Structur 3 bridge in situ

Clip

 Light-curing filling material for temporary restorations


Source: Prof. Dr. Jürgen Manhart, Munich / Germany

- Easy to place easy to remove in one piece
- No damage to preparation boundaries
- Time-saving through light-curing



REF 1285 Syringe 3×4 g

REF 1287 Caps 25 x 0.2 g



Core build-up made easy

If the natural tooth still has sufficient substance, preservation should usually be the goal. With our core build-up classic Rebilda DC and the matching composite root posts Rebilda Post, we offer a core build-up system that behaves like the natural tooth, even under load.

Rebilda® Post

- Glass fibre-reinforced composite post
- Dentine-like elasticity behaviour, high transverse strength
- High level of radiopacity (350 %AI)
- Anatomical shape



REF 1782 Set II 5 posts each of (Ø 1.0 mm, Ø 1.2 mm, Ø 1.5 mm, Ø 2.0 mm), 1 drill each of (Ø 1.0 mm, Ø 1.2 mm, Ø 1.5 mm, Ø 2.0 mm), Ceramic Bond bottle 5 ml, Futurabond U SingleDose 20 pcs., Rebilda DC QuickMix syringe 10 g dentine, accessories

Rebilda® DC

- Dual-curing flowable core build-up and post luting system
- Light-curing for time-saving application
- Additional chemical curing for reliable processing
- Available in 10 g QuickMix syringes and 50 g cartridge



REF 1396 Cartridge 50 g dentine, accessories

REF 1397 Cartridge 50 g blue, accessories

REF 1398 Cartridge 50 g white, accessories

REF 1403 QuickMix syringe 10 g dentine, accessories

REF 1404 QuickMix syringe 10 g blue, accessories

REF 1405 QuickMix syringe 10 g white, accessories



THE global success product

Profluorid Varnish is indispensable for tooth desensitisation. The varnish is easy to apply and dries within seconds. It not only feels good, but also tastes good – because with seven attractive flavours there is definitely one for every patient. Also available in tubes and cartridges.

VOCO Profluorid® Varnish

- Fast desensitisation and fluoride release (5 % NaF ≜ 22,600 ppm fluoride)
- Easy to apply moisture tolerant
- Available in seven pleasant flavours



	SingleDose	Tube
Flavour	$200 \times 0,40$ ml	10 ml
Melon	REF 1271	REF 1267
Caramel	REF 2229	REF 2232
Cherry	REF 2230	REF 2233
Mint	REF 2231	REF 2234
Bubble Gum	REF 2239	-
Cola Lime	REF 2241	REF 2236
Pina Colada	REF 2243	REF 2237
Flavour	4 × 10 ml Tube	
Mixed	REF 2971	

Further trading forms and single shades can be found under www.voco.dental or in the main catalogue







The golden standard of fluoridation

Bifluorid 10 is the daily helper when it comes to treating hypersensitivities. The varnish is versatile and easy to apply – also ideal for orthodontic therapies and bleaching treatments.

Bifluorid 10®

- Fluoride varnish for the treatment of dental hypersensitivity
- Immediately desensitising
- Formation of a protective coat against thermal and mechanical influences
- Special varnish base reinforces long-term effect and deep fluoridation



REF 1616 Bottle 10 g

REF 1617 Bottle 3×10 g, Solvent bottle 10 ml, Pele Tim No. 1

REF 1618 SingleDose 50 pcs., accessories





REF 1619 SingleDose 200 pcs., accessories

Strong fillings without bonding

lonolux®

- Light-curing glass ionomer restorative material
- Excellent working time setting time individually adjustable by light-curing
- No need for conditioning of dental hard tissue
- Does not stick to the instrument, easy to model



REF 2115 Set application capsule 50 pcs. $(10 \times A1, 10 \times A2, 30 \times A3)$

REF 2122 Application capsule 150 pcs. A2 REF 2123 Application capsule 150 pcs. A3

Colourful for children – proven for years

Twinky Star®

- Coloured light-curing componer restorative material with glitter effect
- Eight attractive colours with glitter effect
- Easy and quick application from Caps
- Outstanding polishability



REF 1680 Set Caps 40×0.25 g (5 of each shade), shade guide

















Daily helper for your practice

Fissurit® F

- Light-curing fissure sealant with fluoride, white
- Quick and direct application from the non-dripping NDT®-syringe
- White for application control
- Excellent flow properties and low viscosity



REF 1292 Syringe 2 x 2 g, accessories

REF 1293 Set syringe 5×2 g, accessories

Pele Tim

 Foam pellets for application of liquid or paste preparations without pressure

• 5 practical sizes

• Painfree application without pressure


• Economical application



REF 2250 No. 0 3.000 x ø 3 mm

REF 2252 No. 1 3.000 x ø 4 mm

REF 2253 No. 2 $1.000 \times \emptyset$ 5 mm REF 2254 No. 3 $500 \times \emptyset$ 8 mm REF 2255 No. 4 $500 \times 6 \times 8$ mm



Perfect Bleach – naturally white

Effective, gentle and long-lasting: With our home bleaching system Perfect Bleach you can offer your patients convenient and contemporary tooth whitening that meets the highest aesthetic demands – and this at an attractive price-performance ratio.



Perfect Bleach®

- Simple user-friendly application
- Gentle on tooth substance and existing restorations
- Significant whitening effect already after a short time



REF 1660 Set syringe 6 x 2,4 ml with 10 % carbamide peroxide, Block Out Gel LC syringe 1.2 ml, material for 2 trays, storage box, cosmetic bag, shade guide, accessories

REF 1664 Set syringe 6×2.4 ml with 16 % carbamide peroxide, Block Out Gel LC syringe 1.2 ml, material for 2 trays, storage box, cosmetic bag, shade guide, accessories

Further trading forms and single shades can be found under www.voco.dental or in the main catalogue

Lasts – without etching and bonding!

For a secure adhesion to tooth and restoration: With Bifix SE, the dual-curing, self-adhesive luting system, indirect restorations and root posts can be luted fast and reliably. The composite-based material is suitable for ceramics, zirconium dioxide, composite as well as metal. Bonding agents or conditioner are not needed as all necessary components are already included.

Bifix® SE

- · Fast working: no etching, no bonding
- Odourless
- Easy handling

¹Further trading forms and single shades can be found under www.voco.dental or in the main catalogue



REF 1784 Set QuickMix syringe 3 × 5 g (universal, transparent, white-opaque), accessories

QuickMix syringe 5 g, accessories¹

Denture relining - know how!

Perfectly relined – in one session only. Ufi Gel hard C and Ufi Gel SC are the ideal dream team for the daily prosthetic work. The system offers each application case a special relining solution and can be applied easily and time-saving chairside. Should further adjustments be necessary later, these additions can be made over the years.

Ufi Gel® hard C

- · Hard direct denture relining
- Complete, hard relining chairside, in just one session
- Methyl methacrylate-free
- Neutral taste and odour



REF 2215 Set cartridge 80 g, adhesive bottle 10 ml, accessories

Ufi Gel® SC

- Permanently soft, cold-curing relining material on silicone base
- Top quality relining silicone
- Complete relining in one single session (chairside)
- Remains permanently soft



REF 2041 Set cartridge 50 ml, adhesive bottle 10 ml, glazing base 10 ml, glazing catalyst 10 ml, accessories











 Available in three abrasion grades, thus can be individually adjusted to each situation

Flavour	Coarse (red) Tube 100 g	Medium (yellow) Tube 100 g	Fine (green) Tube 100 g		
Mint	REF 2092	REF 2091	REF 2090		
Caramel	REF 6092	REF 6091	REF 6090		
Cherry	REF 6100	REF 6099	REF 6098		



Your order

Offers are only valid at participating retail partners. Special prices offered by your depot are always taken into consideration.

Errors and omissions excepted. No liability for printing errors.

VOCO customer service

Freecall: 00 800 44 444 555 E-Mail: service@voco.de

handed out by:

This brochure was







