


# Aesthetic, precise, individual

Dentures with the CediTEC System

CediTEC® Denture Base



## CediTEC® *Denture Base*

HIGH-IMPACT PMMA FOR THE FABRICATION OF  
DENTURE BASES FOR REMOVABLE DENTURES



## CediTEC® DB Denture Base

IDEAL IN COMBINATION WITH CEDITEC DENTURE TEETH AND CEDITEC ADHESIVE

### Indications

Fabrication of denture bases for removable dentures

### Advantages

- PMMA with high impact resistance and dentures with the CediTEC System strength for a long time wearing period
- High-quality denture bases, reproducible at any time
- Three shades for an aesthetic result
- Very good polishability – individualisations possible at any time
- Can be processed with all milling units that support a PMMA milling programme

### Presentation

- REF 6195 Disc pink 30 mm, ø 98 mm  
 REF 6196 Disc dark-pink 30 mm, ø 98 mm  
 REF 6197 Disc orange-pink 30 mm, ø 98 mm



### Technical data

Flexural strength*	96 MPa
Modulus of elasticity	2030 MPa
Charpy impact strength + 23 °C	45 kJ/m <sup>2</sup>
Water solubility	< 0,3 µm/mm <sup>3</sup>

\*measured by 3-point flexural strength according to ISO 10477



VOCO GmbH  
 Anton-Flettner-Straße 1-3  
 27472 Cuxhaven  
 Germany  
 Tel.: +49 (0) 4721-719-0  
 Fax: +49 (0) 4721-719-140  
 info@voco.com  
 www.voco.dental

Available from:

