





OUR PARTNERS GRANDIO BLOCS / GRANDIO DISC

	vhf camufacture	www.vhf.de
	Grandio blocs	Grandio disc
Device model	Impression S1 / S2, K4 edition, K5, K5+, R5, Z4, N4	Impression S1 / S2, K4 edition, K5, K5+, R5
Modus	wet + dry	wet + dry
Type of processing	grinding	grinding + milling
Universal holder	yes	yes
Software	DentalCAM 7	DentalCAM 7

	Zirkonzahn	www.zirkonzahn.com
	Grandio blocs	Grandio disc
Device model	Milling Unit M (M1 - M5)	-
Modus	wet	-
Type of processing	grinding	-
Universal holder	yes	no
Software	Zirkonzahn Nesting	-

	imes-icore	www.imes-icore.de
	Grandio blocs	Grandio disc
Device model	CORITEC 245-650i + One	CORITEC 245-650i
Modus	wet	dry
Type of processing	grinding	milling
Universal holder	yes	yes
Software	icam V4.7 / V5 smart	icam V4.7 / V5 smart

	Roland DG	www.rolandeasyshape.com
	Grandio blocs	Grandio disc
Device model	DWX-4, -4W, -51D, -52D, -52DC, -52DCi	DWX-51D, -52D, -52DC, -52DCi
Modus	wet (optional dry)	wet (optional dry)
Type of processing	grinding	grinding + milling
Universal holder	yes	yes
Software	Millbox	Millbox

	exocad GmbH	www.exocad.com
	Grandio blocs	Grandio disc
Device model		
Modus		
Type of processing		
Universal holder		
Software	DentalCAD	DentalCAD

Subject to change without notice. Last update: March 2020