





OUR PARTNERS STRUCTUR CAD

	vhf camufacture	www.vhf.de
	Structur CAD blocs	Structur CAD disc
Device model	Impression S1 / S2, K4 edition, K5, K5+, R5, Z4, N4	Impression S1 / S2, K4 edition, K5, K5+, R5
Modus	wet + dry	wet + dry
Type of processing	milling	grinding + milling
Universal holder	yes	yes
Software	DentalCAM 7	DentalCAM 7

	Zirkonzahn	www.zirkonzahn.com
	Structur CAD blocs	Structur CAD disc
Device model	Milling Unit M (M1 - M5)	–
Modus	wet	–
Type of processing	grinding	–
Universal holder	yes	no
Software	Zirkonzahn Nesting	–

	imes-icore	www.imes-icore.de
	Structur CAD blocs	Structur CAD disc
Device model	CORiTEC 245-650i + One	CORiTEC 245-650i
Modus	wet (category Hybrid Ceramic)	wet + dry (Composite)
Type of processing	grinding	milling
Universal holder	yes	yes
Software	icam V4.7 / V5 smart	icam V4.7 / V5 smart

	Roland DG	www.rolandeasyshape.com
	Structur CAD blocs	Structur CAD disc
Device model	DWX-4, -4W, -51D, -52D, -52DC, -52DCi	DWX-51D, -52D, -52DC, -52DCi
Modus	Composite	Composite
Type of processing	–	–
Universal holder	yes	yes
Software	–	–

Subject to change without notice. Last update: March 2020